
## PRODUCT INFORMATION

**CAUTION:** Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

## CONTENT

Each 1 mL of GELSYN-3 contains:

Sodium Hyaluronate:	8.4 mg
Sodium Chloride:	8.5 mg
Sodium Phosphate, Dibasic:	0.16 mg
Sodium Phosphate, Monobasic:	0.045 mg
Water for Injection:	q.s. to 1.0 mL

## DESCRIPTION

GELSYN-3 is a sterile, buffered solution of highly purified sodium hyaluronate with a molecular weight of approximately 1100 kDa, obtained through fermentation of *Streptococci* of Lancefield group C and is not chemically modified.

## INDICATION

GELSYN-3 is indicated for the treatment of pain in osteoarthritis (OA) of the knee in patients who have failed to respond adequately to conservative non-pharmacologic therapy and simple analgesics (e.g., acetaminophen).

## CONTRAINDICATIONS

- Do not administer to patients with known hypersensitivity (allergy) to sodium hyaluronate preparations.
- Do not inject GELSYN-3 into the knees of patients that have knee joint infections or skin diseases or infections in the area of the injection site.

**Table 1: Primary Endpoint (100 mm WOMAC Pain Subscore)\* for Intent-To-Treat (ITT) Patients**

Variable	GELSYN-3		Commercial Hyaluronan		Difference*	95% CI
	N	Mean ± SE	N	Mean ± SE		
<b>Absolute Values</b>						
Baseline	192	55.2 ± 0.8	188	55.5 ± 0.8	-0.3 ± 1.1	(-2.5, 1.9)
4 Weeks	189	27.0 ± 1.6	183	28.6 ± 1.7	-1.6 ± 1.7	(-4.9, 1.7)
12 Weeks	185	23.6 ± 1.6	178	25.6 ± 1.7	-2.0 ± 1.7	(-5.4, 1.4)
26 Weeks	181	22.2 ± 1.6	175	22.9 ± 1.7	-0.7 ± 1.7	(-4.1, 2.7)
Overall	-	24.3 ± 1.5	-	25.7 ± 1.6	-1.4 ± 1.5	(-4.3, 1.4)
<b>Change from Baseline</b>						
4 Weeks	189	28.1 ± 1.6	183	26.5 ± 1.7	1.6 ± 1.7	(-1.7, 4.9)
12 Weeks	185	31.5 ± 1.6	178	29.5 ± 1.7	2.0 ± 1.7	(-1.4, 5.4)
26 Weeks	181	32.9 ± 1.6	175	32.2 ± 1.7	0.7 ± 1.7	(-2.7, 4.1)
Overall	-	30.8 ± 1.5	-	29.4 ± 1.6	1.4 ± 1.5	(-1.4, 4.3)

\* GELSYN-3 minus Commercial Hyaluronan ± SE. Mixed models include factors for treatment, visit, treatment\* visit interaction, center and baseline level.

## WARNINGS

- Do not concomitantly use disinfectants containing quaternary ammonium salts for skin preparation because sodium hyaluronate can precipitate in their presence.
- Inject into the synovial space only. Do not inject by intravascular route.
- Do not inject outside the synovial space or into the synovial tissue or capsule. An extra-articular injection of the product can cause local adverse events.

## PRECAUTIONS

### General

- The safety and effectiveness of GELSYN-3 in locations other than the knee, and for conditions other than osteoarthritis, have not been established.
- Strict aseptic administration technique must be followed.
- STERILE CONTENTS. EXTERIOR OF SYRINGE IS NOT STERILE. The syringe is intended for single use. The contents of the syringe must be used immediately after its packaging is opened. Do not re-sterilize the product.
- Store in original packaging (protected from light) at room temperature < 25°C (77°F). DO NOT FREEZE.

Do not use GELSYN-3 if package is opened or damaged.

- GELSYN-3 is sensitive to light, and should therefore be used immediately after removal from the carton box.
- Do not use GELSYN-3 in case of severe intra-articular effusion.
- Remove synovial fluid or effusion before each GELSYN-3 injection.
- GELSYN-3 should be used with caution when there is evidence of lymphatic or venous stasis in the leg to be treated.
- The effectiveness of repeat treatment cycles of GELSYN-3 has not been established.

Adverse experience data from clinical use does not show an increased safety risk from retreatment with GELSYN-3. The frequency and severity of adverse events occurring during repeat treatment cycles did not increase over that reported for a single treatment cycle.

## INSTRUCTION FOR PATIENTS

- Provide patients with a copy of the Patient Labeling prior to use.
- Transient pain, sensation of heat, reddening or swelling may occur at the injection site after intra-articular injection of GELSYN-3.
- As with any invasive joint procedure, it is recommended that the patient avoid any strenuous activities or prolonged weight-bearing activities such as jogging or tennis within 48 hours following the intra-articular injection.

## USE IN SPECIFIC POPULATIONS

**Pregnancy:** The safety and effectiveness of GELSYN-3 have not been established in pregnant women.

**Nursing Mothers:** It is not known if GELSYN-3 is excreted in human milk. The safety and effectiveness of GELSYN-3 have not been established in lactating women.

**Children:** The safety and effectiveness of GELSYN-3 have not been established in children.

## ADVERSE REACTIONS

The most common adverse events related to GELSYN-3 injection reported in the clinical study are the following:

- Injection site pain (0.5%)**  
All adverse events related to GELSYN-3 injection reported in the clinical study are provided in the Adverse Events Summary (**Table 4**).
- Post-Marketing Experience**  
Ten adverse events have been reported since GELSYN-3 was first commercialized in Western Europe in 2002, including knee cooling, metallic taste, skin eruption, inflammatory symptoms, burning and itching sensations, tiredness, fever, psoriasis, giant urticaria and edema in both legs, and gouty arthritis. In all the cases, the patients recovered without sequelae.

**Table 2: Secondary Outcome Variables for ITT Patients at 26 Weeks**

Variable	GELSYN-3 (N=192)		Commercial Hyaluronan (N=188)		p
	N	Mean ± SE or %	N	Mean ± SE or %	
1. WOMAC Function	181	28.3 ± 1.7	175	28.0 ± 1.8	0.860
2. WOMAC Stiffness	181	25.9 ± 1.8	175	25.2 ± 1.9	0.735
3. WOMAC Total Score	181	29.0 ± 1.6	175	28.6 ± 1.7	0.811
4. Lequesne Pain	181	2.08 ± 0.16	175	1.79 ± 0.17	0.097
5. Lequesne Walking	181	0.79 ± 0.10	175	0.66 ± 0.11	0.262
6. Lequesne Daily Living	181	1.23 ± 0.12	175	1.10 ± 0.13	0.320
7. Lequesne Total Score	181	4.07 ± 0.32	175	3.53 ± 0.33	0.124
8. Patient Global Pain	181	37.2 ± 2.2	175	33.6 ± 2.3	0.138
9. Patient Global Status	180	25.4 ± 2.0	175	25.7 ± 2.1	0.892
10. Paracetamol Usage	187	46.8 ± 9.7	183	62.8 ± 10.1	0.090
11. Patient Satisfaction	181	82.9%	174	77.0%	0.185
12. Investigator Global Status	181	85.0%	174	76.5%	0.043
13. OMERACT-OARSI	181	89.9%	175	87.7%	0.504

1 to 9: 26-Week Change from Baseline Score (higher=better). Mixed models include factors for treatment, visit, treatment\*visit interaction, center and baseline level.

10: Number of paracetamol rescue medication tablets over the entire 26 Weeks. Mixed model includes factors for treatment and center.

11: Patients satisfied or very satisfied at 26 Weeks. Generalized Estimating Equation (GEE) logistic regression model includes factors for treatment, visit, treatment\*visit interaction and center.

12: Status good or very good at 26 Weeks. Generalized Estimating Equation (GEE) logistic regression model includes factors for treatment, visit, treatment\*visit interaction, center and baseline level.

13: OMERACT-OARSI Success at 26 Weeks. Generalized Estimating Equation (GEE) logistic regression model includes factors for treatment, visit, treatment\*visit interaction and center.

## Potential Adverse Events

The following adverse events are among those that may occur in association with intra-articular injections:

- Arthralgia
- Joint effusion
- Joint warmth
- Arthritis
- Gait disturbance
- Joint stiffness
- Joint swelling
- Injection site pain
- Arthropathy

According to post-marketing experience of other sodium hyaluronate preparations, anaphylactic / anaphylactoid reactions accompanied by transient hypotension (sudden drop in blood pressure), have been rarely reported worldwide, all of which resolved either spontaneously or after conservative treatment.

## Pivotal Clinical Trial

**Study Design:** The safety and efficacy of GELSYN-3 was assessed in a prospective, randomized, double-blind, active control (commercial hyaluronan), non-inferiority study conducted at 23 centers in Europe (Czech Republic, France, Italy, Switzerland, Slovakia, and Germany). A total of 380 patients were enrolled. Patients were given 2 mL intra-articular injections of the randomly assigned GELSYN-3 or commercial hyaluronan once a week for three consecutive weeks, with follow-up visits scheduled for weeks 4, 12 and 26. The primary efficacy variable for this study was the Western Ontario McMaster Universities (WOMAC) pain subscore at week 26 which was required to meet a delta of 8 mm, with secondary efficacy variables including WOMAC total score and stiffness and function subscores, Lequesne Algofunctional Index, global pain assessed by patient, global status assessed by patient and Investigator, paracetamol consumption, patient satisfaction, and overall clinical response based on OMERACT-OARSI criteria. Safety variables included adverse events, pain and local tolerability at the injection site, and global tolerability as assessed by both patient and Investigator.

**Patient Population:** The average age of the 380 intent-to-treat (ITT) patients, defined as receiving at least one intra-articular injection, was 65.0 years, the majority were female (72.9%), BMI ranged from 19.8 to 34.7 kg/m<sup>2</sup>, and the mean duration of osteoarthritis in the target knee was 7.61 years. For 52.6% of the patients, the right knee was the target knee, but 66.1% also had osteoarthritis in the contralateral knee. The anatomical location of target knee osteoarthritis was usually in the medial tibio-femoral region (82.9%), with some degree of joint space narrowing being reported in nearly all patients (95.0%). Patient assessed global pain scores at the screening and baseline visits averaged 65.3 and 65.6, respectively, while the other indices of disease severity at these time-points (i.e. global status, range of motion, WOMAC & Lequesne scoring) were also suggestive of mild to moderate target knee osteoarthritis. The only baseline variable reaching statistical significance with regard to differences between the two treatment groups was venous insufficiency which was more prevalent for commercial hyaluronan patients (5.3% vs 1.0%, p=0.019), with hypercholesterolemia and the requirement for ambulatory assistance also being somewhat (but not statistically significantly) greater (11.7% vs 6.3%, p = 0.073 and 5.9% vs 2.1%, p = 0.069, respectively).

**Efficacy Data:** For the primary outcome measure (change from baseline), the protocol-defined 8 mm non-inferiority margin was met for all time points. In addition, the 95% lower-bound confidence interval of the difference (GELSYN-3 minus commercial hyaluronan) in pain subscore reduction from baseline for the overall 26 week WOMAC pain subscore for the ITT patient population was -1.4 (see Change from Baseline – Overall, **Table 1**). Overall WOMAC pain subscore mean reduction from baseline was 30.8 mm (56%) for the GELSYN-3 treatment group, in contrast to 29.4 mm (53%) for patients receiving commercial hyaluronan.

**Table 3: ITT Patients – Adverse Events\***

Variable	GELSYN-3 (N=192)			Commercial Hyaluronan (N=188)			p Values†
	n*	n	Percent	n	n	Percent	
Any Adverse Event	222	83	43.2%	220	77	41.0%	0.679
<b>Gastrointestinal Disorders</b>							
Toothache	8	7	3.6%	3	3	1.6%	0.337
Diarrhea	1	1	0.5%	4	3	1.6%	0.368
Abdominal pain upper	0	0	0.0%	2	2	1.1%	0.244
<b>General Disorders/Administrative Site Conditions</b>							
Pyrexia	6	5	2.6%	7	5	2.7%	1.000
Injection site pain	3	1	0.5%	4	3	1.6%	0.368
<b>Hepatobiliary disorders</b>							
Cholecystitis acute	0	0	0.0%	2	2	1.1%	0.244
<b>Infections and infestations</b>							
Nasopharyngitis	27	19	9.9%	22	15	8.0%	0.591
Respiratory tract infection	4	4	2.1%	3	3	1.6%	1.000
Urinary tract infection	4	4	2.1%	3	3	1.6%	1.000
Influenza	6	3	1.6%	3	3	1.6%	1.000
Cystitis	1	1	0.5%	5	4	2.1%	0.211
Herpes zoster	0	0	0.0%	2	2	1.1%	0.244
Viral infection	1	1	0.5%	2	2	1.1%	0.620
<b>Musculoskeletal/Connective Tissue Disorders</b>							
Back pain	25	19	9.9%	39	26	13.8%	0.268
Arthralgia	29	18	9.4%	28	22	11.7%	0.506
Pain in extremity	9	9	4.7%	4	2	1.1%	0.062
Musculoskeletal pain	6	5	2.6%	2	2	1.1%	0.449
Neck pain	3	1	0.5%	5	4	2.1%	0.211
Joint swelling	2	2	1.0%	5	2	1.1%	1.000
<b>Nervous System Disorders</b>							
Headache	29	13	6.8%	31	18	9.6%	0.353
<b>Respiratory/Thoracic/Mediastinal Disorders</b>							
Oropharyngeal pain	3	3	1.6%	2	2	1.1%	1.000

\*Adverse events are included in this table only if the incidence in one of the treatment groups exceeded 1%.  
†First n is number of events, and second n is number of patients experiencing an event.  
‡Fisher's exact test.

For the secondary study variables, because this was a non-inferiority study not designed to test superiority, no claims can be made about the statistical significance of any intergroup differences such as those observed at 26 weeks following either three GELSYN-3 or commercial hyaluronan injections (**Table 2**).

**Safety Data: Adverse Events –** Of the 380 patients in the intent-to-treat (ITT) patient population, one or more adverse events were recorded for 160 (42.1%) sometime over the course of the study following the first injection of the assigned hyaluronic acid preparation, by far the most common being back pain (11.8%), arthralgia (10.5%), nasopharyngitis (8.9%), and headache (8.2%). Back pain, arthralgia, and headache were more common in the commercial hyaluronan treatment group than in the GELSYN-3 treatment group (**Table 3**). Adverse events judged to be related to treatment, severe and/or serious were relatively rare (1.6%), and none of the serious adverse events were thought to be treatment-related (**Table 4**). Overall adverse event rates for the GELSYN-3 treatment group were comparable to those of the commercial hyaluronan treatment group.

**Pain at Injection Site/Local Tolerability –** Pain associated with the initial injection of the experimental product or comparator in ITT patients averaged 2.95 on a 10 point scale, declining slightly to 2.80 then 2.65 after the second and third injections, respectively, with local tolerability being judged as good/very good by 91% to 93% of the patients at each of the visits following product administration.

**Patient/Investigator Global Tolerability –** Patient and investigator assessed global tolerability at the 4, 12 and 26-week follow-up visits were nearly identical and comparable to local tolerability scoring, as good/very good scores were obtained from 93% to 97% of the respondents. There were no statistically significant intergroup differences for any of these variables.

**Table 4: ITT Patients – Adverse Events Based on Relationship to Treatment, Severity, Seriousness**

Variable	GELSYN-3 (N=192)			Commercial Hyaluronan (N=188)			p Values†
	n*	n	Percent	n	n	Percent	
<b>Adverse Events Certainly or Probably Related to Treatment</b>							
Any Adverse Event	1	1	0.5%	5	4	2.1%	0.211
<b>General disorders and administration site conditions</b>							
Injection site hematoma	0	0	0.0%	1	1	0.5%	0.495
Injection site pain	1	1	0.5%	1	1	0.5%	1.000
<b>Musculoskeletal and connective tissue disorders</b>							
Arthralgia	0	0	0.0%	1	1	0.5%	0.495
Joint swelling	0	0	0.0%	2	1	0.5%	0.495
<b>Severe Adverse Events</b>							
Any Adverse Event	1	1	0.5%	6	6	3.2%	0.065
<b>Gastrointestinal Disorders</b>							
Diverticulum intestinal	0	0	0.0%	1	1	0.5%	0.495
<b>Hepatobiliary disorders</b>							
Cholecystitis acute	0	0	0.0%	1	1	0.5%	0.495
<b>Infections and infestations</b>							
Influenza	0	0	0.0%	1	1	0.5%	0.495
<b>Injury, poisoning and procedural complications</b>							
Radius fracture	0	0	0.0%	1	1	0.5%	0.495
<b>Musculoskeletal and connective tissue disorders</b>							
Arthritis	0	0	0.0%	1	1	0.5%	0.495
Intervertebral disc protrusion	1	1	0.5%	0	0	0.0%	1.000
<b>Surgical and medical procedures</b>							
Hip arthroplasty	0	0	0.0%	1	1	0.5%	0.495
<b>Serious Adverse Events</b>							
Any Adverse Event	2	2	1.0%	5	4	2.1%	0.445
<b>Gastrointestinal disorders</b>							
Abdominal wall hematoma	1	1	0.5%	0	0	0.0%	1.000
Diverticulum intestinal	0	0	0.0%	1	1	0.5%	0.495
<b>General disorders and administration site conditions</b>							
Pyrexia	0	0	0.0%	1	1	0.5%	0.495
<b>Hepatobiliary disorders</b>							
Cholecystitis acute	0	0	0.0%	1	1	0.5%	0.495
<b>Musculoskeletal and connective tissue disorders</b>							
Arthritis	0	0	0.0%	1	1	0.5%	0.495
Intervertebral disc protrusion	1	1	0.5%	0	0	0.0%	1.000
<b>Surgical and medical procedures</b>							
Hip arthroplasty	0	0	0.0%	1	1	0.5%	0.495

\* First n is number of events, and second n is number of patients experiencing an event.

† Fisher's exact test.

Symbols on packaging	
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.
	Temperature limitation.
	Caution, see insert.
	Protect contents of syringe from sunlight.
	Do not reuse
	Do not use if package damaged.
	Expiration Date
	Catalog number
	Lot number
	Manufacturer
	Sterile – the contents of the syringe have been sterilized by using moist heat. The exterior of the syringe is not sterile.
	Sterile – needle sterilized by ethylene oxide.
	Do not resterilize.

## DETAILED DEVICE DESCRIPTION

GELSYN-3 consists of a buffered physiological solution of 0.84% sodium hyaluronate (**Table 5**).

**Table 5: GELSYN-3 Components (per syringe)**

Component Name	Unitary Amount
Sodium Hyaluronate	16.8 mg
Sodium Chloride	17.0 mg
Sodium Phosphate Dibasic Anhydrous	0.32 mg
Sodium Phosphate Monobasic Dehydrate	0.090 mg
Water for Injection	q.s to 2 mL

## INTERACTIONS

None currently known.

## HOW SUPPLIED

GELSYN-3 is supplied as a sterile solution in a pre-filled glass

## PATIENT INFORMATION

### GELSYN-3 (sodium hyaluronate 0.84%)

Be sure to read the following important information carefully. This information does not take the place of your doctor’s advice. Your doctor has determined that the knee pain you are experiencing is caused by osteoarthritis and that you are a candidate for a non-surgical, non-pharmacological, pain-relieving therapy called GELSYN-3. If you do not understand this information or want to know more, ask your doctor.

## Glossary of Terms

**Hyaluronate:** Hyaluronate is a natural substance found in the human body and is present in very high amounts in joints. The body’s own hyaluronate acts like a lubricant and shock absorber in the joint and is needed for the joint to work properly.

**Non-steroidal anti-inflammatory drug:** Non-steroidal anti-inflammatory drugs are often abbreviated to “NSAIDs”. NSAIDs are drugs, such as aspirin and ibuprofen, for reducing pain, fever and inflammation.

**Osteoarthritis (OA):** Osteoarthritis is a condition that involves the wearing down of cartilage (the protective covering on the ends of your bones) and loss of cushioning fluid in the joint.

## WHAT IS GELSYN-3?

GELSYN-3 contains sodium hyaluronate (0.84% sodium hyaluronate), a viscoelastic material made from bacterial fermentation and is the same material as a natural substance found throughout the body, including joints where it functions as lubricant and shock absorber. Osteoarthritis (pronounced os-TE-o-ar-THRI-tis) (OA) is a type of arthritis that involves the wearing down of cartilage, the protective covering on the end of your bones. In OA, there may not be enough sodium hyaluronate, and there may be a decrease in the quality of the sodium hyaluronate in the joint.

GELSYN-3 is provided in a 2.25 mL glass syringe (half a teaspoon) containing 2.1 mL of solution. The content of the syringe is sterile and is injected directly into the knee, with 2.0 mL being delivered.

## WHAT IS GELSYN-3 USED FOR?

GELSYN-3 is used to relieve knee pain due to OA. It is given to patients who do not get enough relief from non-steroidal anti-inflammatory drugs (NSAIDs) or from simple pain medications, such as acetaminophen, or from exercise and physical therapy.

## HOW IS GELSYN-3 GIVEN?

Your doctor, or other qualified health professional, will inject GELSYN-3 (2 mL) into your knee for a total of three weekly injections.

## ARE THERE ANY REASONS WHY I SHOULD NOT RECEIVE GELSYN-3?

Your doctor will determine if you are a candidate for GELSYN-3 treatment, but you should also be aware that GELSYN-3 should not be administered to patients who:

- have ever had an allergic response to hyaluronate-containing products such as a rash, itching, hives, flushing, swelling of the face, tongue or throat, and/or difficulty breathing;
- have a knee joint infection or skin disease, or infection around the area where the injection will be given, or circulatory problems in the legs.

## WHAT SHOULD MY DOCTOR WARN ME ABOUT?

The following are important treatment considerations for you to discuss with your doctor and understand in order to help avoid unsatisfactory results and complications:

- GELSYN-3 is only for injection into the knee, performed by a qualified doctor.

- GELSYN-3 has not been tested to show better pain relief or safety when combined with other injected medicines.

- Tell your doctor if you are allergic to hyaluronate products.

- For 48 hours after you receive the injection, you should avoid any strenuous activities (such as jogging, tennis, other active sports, heavy lifting) and prolonged weight-bearing activities such as standing on your feet for more than one hour.

- The effectiveness of repeat treatment cycles of GELSYN-3 has not been established.

- Adverse experience data from clinical use does not show an increased safety risk from retreatment with GELSYN-3. The frequency and severity of adverse events occurring during repeat treatment cycles did not increase over that reported for a single treatment cycle.

- Use of GELSYN-3 in joints other than the knee and for conditions other than OA has not been tested.

- GELSYN-3 has not been tested in pregnant or nursing women. You should tell your doctor if you think you are pregnant or if you are nursing a child.

- GELSYN-3 has not been tested in children (≤ 21 years of age).

## WHAT ARE THE POSSIBLE SIDE EFFECTS?

- Extra-articular seepage of GELSYN-3 may cause undesired effects locally.

- Some side effects (also called reactions) may occur during the use of GELSYN-3, with symptoms such as knee pain, stiffness, effusion, or swelling, the sensation of heat, reddening, swelling appearing at the injection site, arthritis, or gait disturbance. These secondary reactions can be relieved by applying ice to the treated joint.

- If any of these symptoms or signs appear after you are given GELSYN-3 or if you have any other problems, you should call your doctor.

## WHAT ARE THE POTENTIAL BENEFITS OF GELSYN-3?

A clinical study involving 380 patients with knee pain due to OA was performed at 23 centers in Europe (Czech Republic, France, Italy, Switzerland, Slovakia, and Germany). The study investigated the safety and effectiveness of GELSYN-3. Patients with osteoarthritic knee joint pain, who had not obtained pain relief with other medications, received either three injections of GELSYN-3 or a commercial hyaluronate into the knee joint. Pain of the knee joint was measured at various times over 4, 12 and 26 weeks. The patients given three injections of GELSYN-3 had the same pain relief as the patients given three injections of commercial hyaluronate for up to 26 weeks after the first injection.

## WHAT DID CLINICAL STUDY SHOW?

The clinical study was conducted at 23 centers in Europe (Czech Republic, France, Italy, Switzerland, Slovakia, and Germany). A total of 380 patients with osteoarthritis were given 2 mL intra-articular injections of either GELSYN-3 or commercial hyaluronate, once a week for three consecutive weeks, with follow-up visits scheduled for weeks 4, 12, and 26. A total of 380 patients were treated, with 192 being injected three times with GELSYN-3 and the remaining 188 being injected three times with a commercially available hyaluronan. Patients were asked to rate their pain under the following five conditions: walking on a flat surface, walking up or down stairs, at night in bed, sitting or lying down (at rest), and while standing. Patients rated their pain from 0 (no pain) to 100 (bad pain) by marking on a 100 mm line. Pain was evaluated in this manner at 1, 4, 12 and 26 weeks after the three injections. The pain scores were used to compare the effectiveness of the GELSYN-3 injections to commercial hyaluronate, and demonstrated that the patients receiving GELSYN-3 experienced the same amount of improvement in knee pain over 26 weeks as those who received the commercially available hyaluronan. Overall

pain score mean reduction from baseline was 30.8 mm (56%) for the GELSYN-3 treatment group, whereas that for the commercial hyaluronan group was 29.4 mm (53%).

## WHAT ADVERSE EVENTS WERE OBSERVED IN THE CLINICAL STUDY?

Of the 380 patients in the study, one or more adverse events were recorded for 160 (42.1%) sometime over the course of the study following the first injection of the assigned hyaluronic acid preparation, by far the most common being back pain (11.8%), arthralgia (10.5%), nasopharyngitis (8.9%), and headache (8.2%). Back pain, arthralgia, and headache were more common in the commercial hyaluronan treatment group than in the GELSYN-3 treatment group. Adverse events judged to be related to treatment, severe and/or serious were relatively rare (1.6%), and none of the serious adverse events were thought to be treatment-related. Overall adverse event rates for GELSYN-3 treated patients were comparable to those of the commercially available hyaluronan group.

## WHAT OTHER TREATMENTS ARE AVAILABLE FOR OA?

If you have OA, there are other things you can do besides getting GELSYN-3. These include:

### Non-drug treatments

- Avoiding activities that cause knee pain

- Exercise

- Physical therapy

- Removal of excess fluid from your knee

### Drug therapy

- Pain relievers such as acetaminophen and narcotics

- Drugs that reduce inflammation (signs of inflammation are swelling, pain or redness), such as aspirin and other nonsteroidal anti-inflammatory drugs (NSAIDs) such as ibuprofen and naproxen

- Steroids injected directly into your knee

## THINGS YOU SHOULD KNOW ABOUT GELSYN-3

- GELSYN-3 is only for injection into the knee, performed by a doctor or other qualified health care professional.

- After you receive the injection, you may need to avoid activities such as jogging, tennis, heavy lifting, or standing for a long time for approximately 48 hours.

- If any of the above symptoms or signs appear after you are given GELSYN-3, or if you have any other problems, you should call your doctor.

## HOW DO I GET MORE INFORMATION ABOUT GELSYN-3?

If you have any questions or would like to find out more about GELSYN-3, you may call 1-800-836-4080 (toll free) or 1-919-474-6700.

<b>DISTRIBUTED BY:</b>
Bioventus LLC 4721 Emperor Blvd., Suite 100 Durham, NC 27703 USA
<b>MANUFACTURED BY:</b>
IBSA Farmaceutici Italia S.r.l. Via Martiri di Cefalonia 2 26900 Lodi Italy

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