90-19681-04

DUROLANE® SJ (small joints)

INSTRUCTIONS FOR USE

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DUROLANE SJ (small joints) est intentionné à être utilisé pour l'injection intra-articulaire pour la symptomatic treatment of mild to moderate osteoarthritis that is diagnosed by synovial joint and for following arthroscopic procedures in the treatment of moderate osteoarthritis. It should be injected by an authorized physician, or in accordance with local legislation.

DUROLANE SJ contient 20 mg/ml. of stabilized non-animal hyaluronic acid in buffered physiological sodium chloride solution pH 7. DUROLANE SJ is a sterile, transparent viscoelastic gel supplied in a 1 mL glass vial. The product is sterilized by gamma radiation.

Hyaluronic acid is identical in all living organisms. It is a natural polysaccharide that is present throughout the tissues of the body, with particularly high concentrations in the synovial fluid and the skin. DUROLANE SJ is composed of biosynthetically produced hyaluronic acid which is purified and stabilized. DUROLANE SJ is degraded in the body by the same metabolic pathway as endogenous hyaluronic acid.

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None known.

Warnings

* DUROLANE SJ should not be injected if the synovial joint is infected or severely inflamed.
* DUROLANE SJ should not be injected if there is an active skin disease or infection present at or near the injection site.

* DUROLANE SJ should not be injected intravascularly or parenterally or in the synovial tissue or capsule.

* Do not re-inject DUROLANE SJ as this may damage the product.

Precautions

* DUROLANE SJ should be used with caution in patients with venous or lymphatic stasis present in the legs.
* DUROLANE SJ has not been tested in pregnant or lactating women or in children.

* A separate syringe of DUROLANE SJ must be used for each individual joint to be treated.

* As with any invasive joint procedure there is a small risk of infection.

* DUROLANE SJ should not be injected if the patient is known to be sensitive to hyaluronic acid based products.

* Local anaesthetics should not be used if the patient is known to be allergic or sensitive to local anaesthetics.

* Injecting under fluoroscopic control and with use of a contrast medium should not be made if the patient is known to be allergic or sensitive to the contrast medium.

* In clinical studies, reinjections have not been studied with a short interval between first and second injection.

* Increase in injection pressure may indicate incorrect extra-articular placement of the needle or overfilling of the joint.

* The effectiveness of DUROLANE SJ following arthroscopic procedures for diagnosis or examination purposes only or in absence of concomitant osteoarthritis of the joint has not been established.

* DUROLANE SJ should be used with caution in patients with pre-existing chondrocalcinosis, as injection may lead to an acute attack of the condition.

Adverse Events

The majority of the reported adverse reactions in clinical studies were described as transient pain, swelling and/or stiffness localized to the joint. These adverse reactions were of mild or moderate intensity and only occasionally required treatment with painkillers or NSAIDs.

The use of other hyaluronic acid preparations in other joints did not reveal any additional unique adverse events.

None of the other adverse reactions that have been reported were interpreted as acute inflammatory arthritis or allergic reactions and they did not need medical attention in the form of surgical intervention, systemic or intra-articular steroids or analgesics.

Adverse events must be reported to the local Biventus representative.

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• Guidance of synovial joints is at the discretion of the treating physician.

• Intra-articular disconfort can be minimized by use of topical freezing agents or subcutaneously delivered local anesthetics.

• Image guidance injection should only be performed by physicians experienced in this type of administration.

Additional information for treatment of synovial joints requiring image guidance

• Guidance of synovial joints is at the discretion of the treating physician.

• Intra-articular disconfort can be minimized by use of topical freezing agents or subcutaneously delivered local anesthetics.

• Image guidance injection should only be performed by physicians experienced in this type of administration.

Additional information for treatment postarthroscopy

• Following the arthroscopic procedure, intra-articular injection should be performed outside the sterile field or the exterior of the syringe is not sterile.

• Smaller joints that typically undergo arthroscopic procedures are the elbow, ankle, and wrist joints.

Please inform your patient that:

* Any other invasive joint procedure it is recommended to avoid strenuous activity the first two days after the injection.

* Some transient reactions related to the injection of DUROLANE SJ, such as pain and/or stiffness, may occur during the first week following the injection can be anticipated. If the symptoms last for more than a week a physician should be contacted.

Performance

* Clinical studies of other hyaluronic acid preparations are similar to DUROLANE SJ and the treatment of osteoarthritis with hyaluronic acid is similar to baseline.

* Some transient reactions related to the injection of DUROLANE SJ, such as pain and/or stiffness, may occur during the first week following the injection can be anticipated. If the symptoms last for more than a week a physician should be contacted.

* The half life of the product is approximately four (4) weeks.

How Supplied

DUROLANE SJ is supplied in a 1 mL glass syringe with a Luer-Lok fitting, packed in a blister pack. The contents of the syringe are sterile. The exterior of the syringe is not sterile.

* The syringe and unused material must be discarded immediately after use. Do not reuse.

* The syringe and unused material must be discarded due to risk of contamination of the unused material and the associated risks including infections.

Disposal should be in accordance with accepted medical practice and applicable national, local or institutional guidelines.

Shelf life and Storage

DUROLANE SJ is supplied in a 1 mL glass syringe with a Luer-Lok fitting, packed in a blister pack. The contents of the syringe are sterile. The exterior of the syringe is not sterile.

* The syringe and unused material must be discarded immediately after use. Do not reuse.

* The syringe and unused material must be discarded due to risk of contamination of the unused material and the associated risks including infections.

Disposal should be in accordance with accepted medical practice and applicable national, local or institutional guidelines.

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DUROLANE, Biventus, and the Biventus logo are registered trademarks of Biventus LLC.