DUROLANE® SJ (small joints)

INSTRUCTIONS FOR USE

Contenu chaque ml. contient:

Hyaluronic acid, stabilized 20 mg
Acid hyaluronique, stabilisé 20 mg
Sodium chloride solution, pH 7 q.s.

Description Description

DUROLANE® SJ (small joints) is intended to be used for intra-articular injection for the symptomatic treatment of mild to moderate osteoarthritis that is indicated by synovial joint pain and/or following arthroscopic procedures in case of degenerative osteoarthritis. It should be injected by an authorized physician, or in accordance with local legislation.

DUROLANE SJ contains 20 mg/ml. of stabilized non-animal hyaluronic acid buffered physiological sodium chloride solution pH 7. DUROLANE SJ is a sterile, transparent viscous gel supplied in a 1 mL glass vial.

Hyaluronic acid is identical in all living organisms. It is a natural polysaccharide that is present throughout the tissues of the body, with particularly high concentrations in the synovial fluid and the skin. DUROLANE SJ is composed of biosynthetically produced hyaluronic acid which has been purified and stabilized. DUROLANE SJ is degraded in the body by the same metabolic pathway as endogenous hyaluronic acid.

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DUROLANE SJ should not be injected intravascularly or into the synovial tissue or capsule.

• Do not re-inject DUROLANE SJ as this may damage the product.

Prae-contraindications

DUROLANE SJ should be used with caution in patients with venous or lymphatic stasis present in the legs.

DUROLANE SJ has not been tested in pregnant or lactating women or in children.

• A separate syringe of DUROLANE SJ must be used for each individual joint to be treated.

• As with any invasive joint procedure there is a small risk of infection.

• DUROLANE SJ should not be injected if the patient is known to be sensitive to hyaluronic acid based products.

• Local anaesthetics should not be used if the patient is known to be allergic or sensitive to local anaesthetics.

• Injecting under fluoroscopic control and without use of a contrast medium should not be made if the patient is known to be allergic or sensitive to the contrast medium.

• In clinical studies, rejections have not been studied with a short interval between first and second injections.

• Increase in injection pressure may indicate incorrect extra-articular placement of the needle or overfilling of the joint.

• The effectiveness of DUROLANE SJ following arthroscopic procedures for diagnosis or examination purposes only or in absence of concomitant osteoarthritis of the joint has not been established.

• DUROLANE SJ should be used with caution in patients with pre-existing chondrocalcinosis, as injection may lead to an acute attack of the condition.

Adverse Events

The majority of the reported adverse reactions in clinical studies were described as transient pain, swelling and/or stiffness localized to the joint. These adverse reactions were of mild or moderate intensity and only occasionally required treatment with painkillers or NSAIDs.

The use of other hyaluronic acid preparations in other joints did not reveal any additional unique adverse events.

None of the other adverse reactions that have been reported were interpreted as acute inflammatory arthritis or allergic reactions and they did not need medical attention in the form of surgical intervention, systemic or intra-articular steroids or analgesics.

Adverse events must be reported to the local Bioventus representative.

The safety and effectiveness of DUROLANE SJ concomitantly with other intra-articular injectables have not been established.

Administration

General administration information

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• DUROLANE SJ should be injected using strict aseptic technique.

• DUROLANE SJ should be injected into the joint cavity only.

• Intra-articular injection in certain synovial joints will require image guidance to ensure accurate placement and depth of the needle and orientation of the joint.

• The route for intra-articular injection with or without image guidance should be chosen so that damage to adjacent vital structures is avoided.

• The injection site should be washed with alcohol or other suitable antisepsis before injection.

• Remove joint effusion, if present, before injecting DUROLANE SJ. The same needle should be used for both removal of effusion and injection of DUROLANE SJ.

• The recommended needle size is 18 to 25 G and with adequate length.

• Use of smaller diameter needles increases pressure required to deliver the product.

Additional information for treatment of synovial joints requiring image guidance

• Guidance of synovial joints is at the discretion of the treating physician.

• Injecting discomfort can be minimized by use of topical freezing agents or subcutaneously delivered local anesthetics.

• Image guided injection should only be performed by physicians experienced in this type of administration.

Additional information for treatment post arthroscopy

• Following the arthroscopic procedure, intra-articular injection should be performed outside the sterile field at the exterior of the syringe is not sterile.

• Smaller joints that typically undergo arthroscopic procedures are the elbow, ankle, and wrist joints.

Please inform your patient that:

• All joint therapy joint procedure it is recommended to avoid strenuous activity the first two days after the injection.

• Some transient reactions related to the injection of DUROLANE SJ, such as pain and/or swelling of the joints, may appear within the first week following the injection can be anticipated. If the symptoms last for more than a week a physician should be contacted.

Performance

• Clinical studies of other hyaluronic acid preparations show that DUROLANE SJ is the treatment of osteoarthritis pain, the pain relief is better than baseline over baseline values. Some studies also showed improvements favoring the hyaluronic acid treated group over that of the control therapy, such as saline and corticosteroids. Improvement in pain and functional outcome to 6 months post-treatment were observed.

• The half-life of the product is approximately four (4) weeks.

How Supplied

DUROLANE SJ is supplied in a 1 mL glass vial with a Luer-Lok fitting, packed in a blister pack. The contents of the syringe are sterile. The exterior of the syringe is not sterile.

• The syringe and any unused material must be discarded immediately after use. The exterior of the syringe is not sterile due to risk of contamination of the unused material and the associated risks including infections.

Disposal should be in accordance with accepted medical practice and applicable national, local or institutional guidelines.

Shielding and Storage

DUROLANE SJ should be stored, in its original packaging, up to 30°C. The expiry date is indicated on the package and should not be used beyond that date. Protect from freezing.

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DUROLANE, Biventus, and the Biventus logo are registered trademarks of Biventus LLC.